

Automatic Image Segmentation of Light-Sheet Zebrafish Scans

MICROBS Lab

Supervising Professor :
Selman Sakar

Supervising TA :
Artur Krzysztof Banach


Aly Elbindary

- Background & Objectives
- Methodology
- Results
- Discussion

Background & Objectives


Goals & Motivation

- What do we want to achieve? Why?




Goals & Motivation

- What do we want to achieve? Why?




Methodology

■ Manual Segmentation



Methodology


- Automatic Segmentation – Introduction to SAM



Segment Anything Model (SAM) – Architecture Overview [5]

Methodology

- Automatic Segmentation – Introduction to SAM




Example images with overlaid masks generated by SAM [5]


Methodology

- Raw Data Adaptation

Raw Zebrafish Frame Loaded
Without Adaptation




Raw Mask Frame Loaded
Without Adaptation




Methodology

- Raw Data Adaptation


Before Pixel Conversion




After Pixel Conversion




Raw Mask



90° Rotation Around Center



180° Rotation Around Middle Horizontal Axis



- Loss Function & Evaluation Metrics


Loss Function : Focal Loss [3]

$$CE(p, y) = \begin{cases} -\log(p) & \text{if } y = 1 \\ -\log(1 - p) & \text{otherwise} \end{cases}$$

$$p_t = \begin{cases} p & \text{if } y = 1 \\ 1 - p & \text{otherwise} \end{cases}$$

$$FL(p_t) = (1 - p_t)^\gamma CE(p_t)$$

Evaluation Metric : Intersection over Union (IoU) [7]



$$IoU = \frac{\text{Area of Overlap}}{\text{Area of Union}}$$

Results


Results

- Finetuning SAM – Version : facebook/sam-vit-base (93.7M)


Train Loss




Train IoU



Test Loss




Test IoU



Results


- Finetuning SAM – Version : facebook/sam-vit-base (93.7M)




Results

- Finetuning SAM – Version : facebook/sam-vit-huge (641M)


Train Loss




Train IoU



Test Loss





Test IoU



Results


- Finetuning SAM – Version : facebook/sam-vit-huge (641M)




Results

- Finetuning SAM – Version : wanglab/medsam-vit-base


Train Loss




Train IoU



Test Loss




Test IoU




Results

- Finetuning SAM – Version : wanglab/medsam-vit-base



Results

- Finetuning SAM – Multi-Class Image Segmentation




Discussion

Discussion

■ Results Analysis – Which Model Yields The Better Results?

facebook/sam-vit-huge Model inferences



wanglab/medsam-vit-base Model inferences




■ Drawbacks of SAM & Potential Improvements

Improvements List :

- Hyperparameter Search
- Multi-Class Segmentation Fix
- Data Augmentation & Improvements

Prompt Bounding Box Example



Summary & Conclusion

- [1] Meta AI. *Segment Anything*. 2023.
- [2] Rafael Guedes. *SAM: Segment Anything Model - Quickly customize your product landing page with SAM*. <https://towardsdatascience.com/sam-segment-anything-model-4b25a47245f2>. 2021.
- [3] Shruti Jadon. “A survey of loss functions for semantic segmentation”. In: *2020 IEEE Conference on Computational Intelligence in Bioinformatics and Computational Biology (CIBCB)*. 2020, pp. 1–7. DOI: 10.1109/CIBCB48159.2020.9277638.
- [4] Jan Witowski. *3D Slicer Tutorial*. YouTube playlist. 2019. URL: <https://www.youtube.com/playlist?list=PLeaIM0zU1Eqswa6Pskg9uMq15LiWWYP39>.
- [5] Alexander Kirillov et al. *Segment Anything*. 2023. arXiv: 2304.02643 [cs.CV]. URL: <https://arxiv.org/pdf/2304.02643.pdf>.
- [6] Alexandra L’Heureux et al. “Machine Learning With Big Data: Challenges and Approaches”. In: *IEEE Access* 5 (2017), pp. 7776–7797. DOI: 10.1109/ACCESS.2017.2696365.
- [7] Adrian Rosebrock. *Intersection over Union (IoU) for object detection*. URL: <https://pyimagesearch.com/2016/11/07/intersection-over-union-iou-for-object-detection/>.
- [8] Saskia Dwi Ulfah. *Brain MRI Segmentation with Segment Anything Model (SAM)*. 2023. URL: <https://medium.com/@sdwiulfah/brain-mri-segmentation-with-segment-anything-model-sam-16d0b4101a85>.

Python Functions that Allow Proper Mask Rotation

```
1 """
2 Rotate .nrrd data in different ways
3 """
4 def rotate_around_middle_horizontal_axis(data, angle):
5     # Get the height and width of each frame
6     height, width, num_frames = data.shape
7
8     # Calculate the middle index along the horizontal axis for each
9     # frame
10    middle_horizontal_indices = width // 2
11
12    # Initialize an array to store the rotated data
13    rotated_data = np.empty_like(data)
14
15    # Iterate over each frame
16    for i in range(num_frames):
17        # Rotate the data within the frame around the middle horizontal
18        # axis
19        rotated_frame = np.rot90(data[:, :, i], k=angle // 90, axes=(1,
20        0))
21
22        # Flip the rotated frame along the horizontal axis to align
23        # with the original orientation
24        rotated_frame = np.flip(rotated_frame, axis=1)
25
26        # Store the rotated frame in the output array
27        rotated_data[:, :, i] = rotated_frame
28
29    return rotated_data
30
31 def rotate_data(data, angle):
32     # Rotate the data by the specified angle
33     rotated_data = np.rot90(data, k=angle//90, axes=(0, 1)) # Adjust
34     # axes if needed
35     return rotated_data
36 """
```

Python Function that Allows .tif Pixel Value Conversion

```
1 def adjust_jpg(image):
2     """
3     Takes in a jpeg that has been extracted from a .tif file and
4     adjusts the pixel intensity in order to get a readable image
5
6     Inputs :
7         - image : original jpeg to be adjusted
8     Outputs :
9         - image : adjusted jpeg
10    """
11    # Convert image to uint8 type
12    if image.dtype != np.uint8:
13        image = (image / (np.max(image) / 255)).astype(np.uint8)
14
15    return image
```

Python Code for Multi-Class SAM Model

```
1 class MultiSAM(nn.Module):  
2     def __init__(self, num_classes):  
3         super(MultiSAM, self).__init__()  
4  
5         self.sam = SamModel.from_pretrained("facebook/sam-vit-base")  
6         self.conv_layer = nn.Conv2d(in_channels=1, out_channels=  
7             num_classes, kernel_size=1)  
8  
9         def forward(self, x):  
10             x = self.sam(pixel_values = x["pixel_values"],  
11                 input_boxes = x["input_boxes"],  
12                 multimask_output=False).pred_masks.squeeze(1).to  
13  
14             x = self.conv_layer(x)  
15  
16             return x  
17
```